

Information on Devices Located in the Laboratory			
Device Specifications		Device Supervisor	
Brand / Model	Hemostatic Forceps (straight) (TIPKİMSAN)	Name Surname	Asst. Prof. Necmiye ÇÖMLEKÇİ Res. Asst. Meryem YÜCEL
Fixed Asset Registration Number	253.3.5/22/3316 253.3.5/22/3317	Phone	0 378 501 10 00
Device Status / Quantity	2 piece	E-Mail	ncomlekci@bartin.edu.tr meryemyucel@bartin.edu.tr
Laboratory Where It Is Located	First Aid Laboratory		

Photo of the Device



Technical Specifications:

- It is manufactured using stainless materials.
- The hemostatic forceps are straight.
- It opens and closes smoothly without jamming, and when fully closed, the jaws fit together seamlessly without any gaps.
- There are no sharp edges except for the cutting parts.
- There is absolutely no burr, gap, etc. on the surgical instrument at the joints in the mouth area or on its body.
- It is manufactured to be suitable for all kinds of sterilization conditions.
- There is not the slightest wobble in the fixed joints of the forceps, and the jaw blades and serrations interlock without leaving any gaps.
- The arms of the forceps are flexible like springs.
- Its size is 14 cm $\pm$ 2 cm.