

Information on Devices Located in the Laboratory			
Device Specifications		Device Supervisor	
Brand / Model	Child Basic Life Support Training Manikin with Feedback (LAERDAL)	Name Surname	Asst. Prof. Necmiye ÇÖMLEKÇİ Res. Asst.Meryem YÜCEL
Fixed Asset Registration Number	255.1.5/21/5164	Phone	0 378 501 10 00
Device Status / Quantity	1 piece	E-Mail	ncomlekci@bartin.edu.tr meryemyucel@bartin.edu.tr
Laboratory Where It Is Located	First Aid Laboratory		
Photo of the Device 			
Technical Specifications: - The manikin must have a lower body made of special high-quality plastic, be in the structure of a child half body, and have chest skin and a face with realistic texture. - The manikin must be designed and feature a structure that helps improve both CPR training quality, classroom training efficiency, and students' self-development. - Chest rise can be observed during artificial respiration on the manikin. - For quality CPR, nose pinching, head tilt, chin lift, and jaw thrust can be performed on the manikin. - Correct compressions can be seen and heard. - Anatomically correct oral and nasal passages can be made on the manikin. - Additionally, thanks to its Bluetooth feature, the manikin allows CPR performance evaluation by installing the free software on any phone or tablet PC. - Thanks to this feature, values such as compression rate, compression depth, and ventilation volume during CPR can be observed instantly and live, and can be projected onto a projector or TV in the classroom so everyone can see and to create a CPR competition environment, making training realistic. - Included with the manikin are an upper tracksuit, 2 spare lungs, 2 spare face skins, 10 cleaning wipes, the original user manual, and a soft carrying bag that can also be used as a training mat.			