

Information on Devices Located in the Laboratory			
Device Specifications		Device Supervisor	
Brand / Model	Half-Body Adult Basic Life Support Training Manikin (GANGRENE)	Name Surname	Asst. Prof. Necmiye ÇÖMLEKÇİ Res. Asst.Meryem YÜCEL
Fixed Asset Registration Number	255.1.5.11-111	Phone	0 378 501 10 00
Device Status / Quantity	1 piece	E-Mail	ncomlekci@bartin.edu.tr meryemyucel@bartin.edu.tr
Laboratory Where It Is Located	First Aid Laboratory		

Photo of the Device



Technical Specifications:

- The manikin is produced from special quality plastic in an adult upper body structure and features realistic skin texture on the chest and face.
- The manikin is designed in a structure and with features that help improve both CPR training quality and classroom training efficiency, as well as helping students develop themselves.
- Chest rise can be observed during artificial respiration on the manikin.
- Nose pinching, head tilt, chin lift, and jaw thrust can be performed on the manikin for quality CPR.
- Anatomically correct oral and nasal passages can be made on the manikin.
- The accuracy of compressions can be checked visually and audibly.
- There is a warning system with a light indicating the correct compression depth.