

Information on Devices Located in the Laboratory			
Device Specifications		Device Supervisor	
Brand / Model	Pediatric Basic Life Support Training Manikin (NASCO)	Name Surname	Asst. Prof. Necmiye ÇÖMLEKÇİ Res. Asst.Meryem YÜCEL
Fixed Asset Registration Number	255.8.2/20/4160	Phone	0 378 501 10 00
Device Status / Quantity	1 piece	E-Mail	ncomlekci@bartin.edu.tr meryemyucel@bartin.edu.tr
Laboratory Where It Is Located	First Aid Laboratory		

Photo of the Device



Technical Specifications:

- The manikin has a lower body made of special quality plastic, a child half-body structure, and realistic-textured chest skin and face.
- The manikin is designed and featured to help improve both CPR training quality, classroom training efficiency, and students' self-development.
- Chest rise can be observed during artificial respiration on the manikin.
- On the manikin, nose pinch, head tilt, chin lift, and jaw thrust can be performed for quality CPR.
- The compressions can be seen and heard when performed correctly.
- Anatomically correct oral and nasal passages can be created on the manikin.
- Additionally, thanks to its Bluetooth feature, the manikin allows for CPR performance evaluation by installing the free program on any phone or tablet PC.
- Thanks to this feature, values such as compression rate, compression depth, and ventilation volume during students' CPR can be observed instantly and live, and can be projected onto a projector or TV so that everyone in the classroom can see and to make training realistic by providing a CPR competition environment.
- Along with the manikin, an upper tracksuit, 2 spare lungs, 2 spare face skins, 10 cleaning wipes, the original user manual, and a soft carrying bag that can also be used as a training mat are included.