

Information on Devices Located in the Laboratory			
Device Specifications		Device Supervisor	
Brand / Model	VZN Tweezers (toothed)(TIPKİMSAN)	Name Surname	Asst. Prof. Necmiye ÇÖMLEKÇİ Res. Asst.Meryem YÜCEL
Fixed Asset Registration Number	253.3.5/22/3312	Phone	0 378 501 10 00
Device Status / Quantity	1 piece	E-Mail	ncomlekci@bartin.edu.tr meryemyucel@bartin.edu.tr
Laboratory Where It Is Located	First Aid Laboratory		
Photo of the Device 			
Technical Specifications: • It is manufactured using stainless materials. • The forceps have no play whatsoever at their fixed joints, and their tips and serrations interlock without leaving any gap. • It has serrated jaws. • The gaps between the teeth and serrations are clean. • The arms of the forceps are flexible like springs. • The forceps can be opened and closed easily without sticking, and when fully closed, the jaws of their parts align with each other and seat without leaving any gap. • There are no sharp edges except for the cutting parts. • There is no burr, gap, etc. in any way at the joints and body of the surgical instrument's mouth part. • It is manufactured to be suitable for all kinds of sterilization conditions. • Its size is 14 cm $\pm$ 2 cm.			